# ANTIGENIX AMERICA Inc. 1-800-558-1008 PRODUCT INFORMATION SUMMARY

TMB Substrate Solutions: cat no: ES200

Solutions Part A and Part B 100 mL (2)

## (ELISA Construction Pack- component)

### **DESCRIPTION:**

TMB (3,3', 5,5' -Tetramethyl-Benzidine) substrate solution is used for ELISA applications in immunologic and serologic procedures where HRP (Horseradish Peroxidase)is applied. It is non-carcinogenic.

### HANDLING AND STORAGE:

### **PRODUCT USE:**

Provided as 100 mL Solutions Part A and Part B.

Store DARK at 4 Deg. C.

Prepare substrate solution no more than **15 minutes before last incubation of assay:** Mix one part TMB Solution A with one part TMB Solution B in a clean container. If, upon mixing, TMB soltion turns blue - TMB solution is **contaminated - DO NOT USE**. Use mixed substrate solution WITHIN 2 HOURS, AND AVOID DIRECT LIGHT.

Apply TMB soultion to wells of assay as substrate. For 96-well microplates, usually 100 uL per well is recommended. Be sure the ELISA plate is **thoroughly washed (5 times)** before adding TMB solution. (Wash steps are critical!)

For most ELISA applications, incubate TMB substrate solution for 15-20 minutes at 37 Deg. C. Then add 100 uL stop solution (2N Sulphuric Acid) **and mix by gently tapping plate**, then read O.D. at 450 nm within 30 minutes. - or follow customer's usual protocol.

**CAUTION:** Substrate solution A contains 20% acetone - **flammable -** keep away from sources of heat or flame.