

# Super-X Plex™ Flow Cytometry Assay

## Product Information Summary

**DESCRIPTION:** Flow Cytometry Assay components to measure **human Uteroglobulin (CC10)**, in human Serum, Plasma, Cell Culture Supernatant or Cell Lysate samples. Can be custom Multi-Plexed with other human assays.

Cell Lysate Prep Buffer (**CLP200**) or Sample Dilution Buffer  
( to dilute high concentration samples) - (**SDL200**)  
also available

**Standard:** One vial of lyophilized standard containing  
Pre-mixed multiplexed analytes ( see vial insert) -single  
standard vial can be used for one or more assays indicated

**CAUTION:** Contains Sodium Azide - ( 0.05%) - Dispose contents carefully, flush with water to avoid reaction with lead drain lines.

Recovery 70-130% standard dose

**WWW.ANTIGENIX.com 1-800-558-1008**