

ANTIGENIX AMERICA Inc.

**Super-X Plex™ Flow Cytometry Assay
Multi-Plex Panel (6 Plex)**

Product Information Summary

**Non-Human Primate Th1/Th2 - Cat No.PMX776T
32 tests**

DESCRIPTION: Multi-Plex (6-Plex) Flow Cytometry Assay Panel
- components to measure **Non-human Primate : IL-2, IL-4,
IL-5, IL-6, IFN Gamma, and TNF Alpha simultaneously-** in
Primate Serum, Plasma, Cell Culture Supernatant or Cell
Lysate samples.

NOTE: ORDER a Sample Diluent Kit (to match sample type):

Human/Primate Serum/Plasma (**DKH100**) or
Cell Culture Supernatant (**DKC100**) or
Cell Lysate (**DKL100**)

Cell Lysate Prep Buffer (**CLP200**) or Sample Dilution Buffer
(to dilute high concentration samples) -(**SDL200**)
also available.

(STARTER Kit Included)

CONTENTS:

Antibody Conjugated Beads 1.5 mL (1X) (**Pre-Mixed , 6-Plex**)
Detection Antibody - Biotin 0.45 mL (2X) (**Pre-Mixed , 6-Plex**)
Detection Antibody Diluent 0.45 mL (2X)
Standard: One vial of lyophilized standard containing
Pre-mixed multiplexed analytes (see vial insert) -single
standard vial can be used for one or more assays indicated

STARTER Kit: (Human/Primate Panel)

Reading Buffer 5.0 mL (10X)
Wash Buffer 15.0 mL (10X)
Filter plate with lid (96 wells)
Plate Sealers (6)
Streptavidin-PE 3.0 mL (1X)
PCR 8 tube strip 2 strips

STORAGE: Store Refrigerated (4-8° C. - in dark)

CAUTION: Contains Sodium Azide - (0.05%) - Dispose contents
carefully, flush with water to avoid reaction with lead drain
lines.

PROTOCOL: Refer to Kit protocol provided with assay ordered.
PROTOCOL can also be DOWNLOADED @ www.antigenix.com

Assay Specifications:

Bead Specs:	Sensitivities (pg/mL)	Range: (pg/mL)
S4P3	Primate TNFa < 1.0	2.0 - 5,000
S4P5	Primate IL-5 < 5.0	10.0 - 5,000
S4P7	Primate IL-2 < 5.0	10.0 - 5,000
S4P9	Primate IFN Gamma < 5.0	10.0 - 5,000
S5P4	Primate IL-4 < 1.0	5.0 - 2,000
S5P6	Primate IL-6 < 2.0	5.0 - 5,000

Sample Size: 15.0 uL per test
X-Reactivity within panel : negligible
Intra Assay CV: <10% Inter Assay CV <20%
Recovery 70-130% standard dose

RESEARCH USE ONLY
NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE.

WWW.ANTIGENIX.com 1-800-558-1008