

ANTIGENIX AMERICA Inc.

**Super-X Plex™ Flow Cytometry Assay  
Multi-Plex Panel ( 6 Plex)**

Product Information Summary

**Human Th1/Th2 - Panel #1 Cat No.HMX776  
96 tests**

**DESCRIPTION:** Multi-Plex ( 6 Plex) Flow Cytometry Assay Panel  
- components to measure **human: IL-2, IL-4, IL-6, IL-10,  
IFN Gamma, and TNF Alpha simultaneously-** in human Serum,  
Plasma, Cell Culture Supernatant or Cell Lysate samples.

**NOTE: ORDER a Sample Diluent Kit ( to match sample type):**

Human Serum/Plasma (**DKH100**) or  
Cell Culture Supernatant (**DKC100**) or  
Cell Lysate (**DKL100**)

Cell Lysate Prep Buffer (**CLP200**) or Sample Dilution Buffer  
( to dilute high concentration samples) -(**SDL200**)  
also available.

(**STARTER Kit Included**)

**CONTENTS:**

**Antibody Conjugated Beads** 5.0 mL (1X) (**Pre-Mixed , 6-Plex**)  
**Detection Antibody - Biotin** 1.5 mL (2X) (**Pre-Mixed , 6-Plex**)  
**Detection Antibody Diluent** 1.5 mL (2X)  
**Standard:** One vial of lyophilized standard containing  
**Pre-mixed multiplexed** analytes ( see vial insert) -single  
standard vial can be used for one or more assays indicated

**STARTER Kit: ( Human Panel)**

Reading Buffer 5.0 mL (10X)  
Wash Buffer 15.0 mL (10X)  
Filter plate with lid ( 96 wells)  
Plate Sealers (6)  
Streptavidin-PE 3.0 mL (1X)  
PCR 8 tube strip 2 strips

**STORAGE:** Store Refrigerated ( 4-8° C. - in dark)

**CAUTION:** Contains Sodium Azide - ( 0.05%) - Dispose contents  
carefully, flush with water to avoid reaction with lead drain  
lines.

**PROTOCOL:** Refer to Kit protocol provided with assay ordered.  
PROTOCOL can also be DOWNLOADED @ [www.antigenix.com](http://www.antigenix.com)

**Assay Specifications:**

<b>Bead Specs:</b>	<b>Sensitivities (pg/mL)</b>	<b>Range: (pg/mL)</b>
S4P3	Human TNF Alpha < 1.0	2.0 - 1,000
S4P7	Human IL-2 < 3.0	10.0 - 2,000
S4P9	Human IFN Gamma < 3.0	10.0 - 2,000
S5P4	Human IL-4 < 1.0	3.0 - 1,000
S5P6	Human IL-6 < 5.0	10.0 - 5,000
S5P10	Human IL-10 < 2.0	5.0 - 5,000

Sample Size: 15.0 uL per test  
X-Reactivity within panel : negligible  
Intra Assay CV: <10%      Inter Assay CV <20%  
Recovery 70-130% standard dose

**RESEARCH USE ONLY**  
**NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE.**

[WWW.ANTIGENIX.com](http://WWW.ANTIGENIX.com) 1-800-558-1008