ANTIGENIX AMERICA Inc.

Super-X Plex[™] Flow Cytometry Assay Multi-Plex Panel (5-Plex)

Product Information Summary

Mouse Th1/Th2 Cat No. MMX775 96 tests

DESCRIPTION: Multi-Plex (5-Plex) Flow Cytometry Assay Panel - components to measure Mouse: IL-2, IL-4, IL-5, IFN Gamma, and TNF Alpha simultaneously- in Murine Serum, Plasma, Cell Culture Supernatant or Cell Lysate samples.

NOTE: ORDER a Sample Diluent Kit (to match sample type):

Mouse/Rat Serum/Plasma (DKM100) or Cell Culture Supernatant (DKC100) or Cell Lysate (DKL100)

Cell Lysate Prep Buffer (CLP200) or Sample Dilution Buffer (to dilute high concentration samples) - (SDL200) also available.

(STARTER Kit Included)

CONTENTS:

Antibody Conjugated Beads 5.0 mL (1X) (Pre-Mixed , 5-Plex)
Detection Antibody - Biotin 1.5 mL (2X) (Pre-Mixed , 5-Plex)
Detection Antibody Diluent 1.5 mL (2X)

Standards: One vial of lyophilized standards containing

Pre-mixed multiplexed analytes (see vial insert) -single
standard vial can be used for one or more assays indicated

STARTER Kit: (Mouse Panel)

Reading Buffer 5.0 mL (10X)
Wash Buffer 15.0 mL (10X)
Filter plate with lid (96 wells)
Plate Sealers (6)
Streptavidin-PE 3.0 mL (1X)
PCR 8 tube strip 2 strips

STORAGE: Store Refrigerated (4-8°C. - in dark)

CAUTION: Contains Sodium Azide - (0.05%) - Dispose contents carefully, flush with water to avoid reaction with lead drain lines.

PROTOCOL: Refer to Kit protocol provided with assay ordered.

Assay Specifications:

Bead Class:	Sensitivities (pg	/mL) Range: (pg/mL)
S4P8 Mous S5P2 Mous S5P7 Mous	se IL-2 < 2.0 se IFN Gamma < 0.5 se IL-4 < 0.5 se IL-5 < 5.0 se TNF Alpha < 0.5	5.0 - 5,000 2.0 - 5,000 2.0 - 1,500 10.0 - 5,000 2.0 - 1,500

Sample Size: 15.0 uL per test

X-Reactivity within panel: negligable Intra Assay CV: <10% Inter Assay CV <20% Recovery 70-130% standard dose

RESEARCH USE ONLY NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE.

WWW.ANTIGENIX.com 1-800-558-1008