ANTIGENIX AMERICA Inc.

PRODUCT INFORMATION SUMMARY

Human KLOTHOProduct Number RHF702CKX"SUPER X" ELISA Kit5 Plate Kit (5 x 96 tests)RHF702CKX2 (2 x 96 tests)RHF702CKX1 (single plate)

PRE-COATED ELISA Plates(single, 2 or 5 plates)Biotin-Labeled tracer25.0 ug or 0.5 mL liquid - see vial)Antigen Standard5.0 ugStreptavidin-HRP0.5 mLTMB Substrate25 mL x 2 (15 mL - 2 plate)Wash Buffer Concentrate125 mL (20X) (50 mL - 2 plate)

DESCRIPTION:

Human KLOTHO **SUPER X ELISA Kit** provides single, two or five **PRE-COATED** microplates (coated with antigen affinity purified capture antibody stabilized with our proprietary ELISA coating/blocking reagent. A biotin labeled tracer antibody, antigen standard, HRP developing reagents and wash buffer are included.

Note: Reconstitute components only when ready to run assay.

TRACER ANTIBODY:

Provided as 25 ug (lyophilized)or as 0.5 mL liquid

@ 50.0 ug/mL (see vial) of Biotin labeled, antigen-affinity
purified antibody, additive-free. Reconstitute the 25 ug
(lyophilized vial) in 500 uL sterile water containing 0.1% BSA.

(FREEZE aliquots for long-term storage).

******For liquid vial store refrigerated only (contains preservative) ******.

STANDARD: Provided as 5.0 ug of recombinant Human Klotho. Quickspin and reconstitute in 50 uL of distilled water (pH 7.6) Further dilutions can be made in Diluent provided (ED100) or with 0.1% BSA in PBS.

DEVELOPING REAGENTS:

- Streptavidin-HRP (S100180CX) 0.5 mL store @ -20 Deg. C.
- TMB Substrate Solutions Part A and Part B (25.0 mL each) cat # ES200CX
- Wash Buffer Concentrate (20X concentrate) Mix 1 volume of wash buffer with 19 volumes of distilled water. Stable for one month @ 4 Deg. C. once mixed to working volume.
- Diluent (10X) -Dilute 1 part with 9 parts distilled water
- STOP Solution (1X) **Caution:** Caustic